

SPECTROSCOPY

We will proceed in steps:

1) If a light with power P_d hits the photodiode we will see a voltage

$$V_{pd} = G P_d, \quad G := \text{photodiode gain}$$

2) If a light intensity I_0 enters the cell we will have

$$I_{out} = I_0 e^{-K(\omega) \Delta z} \quad \Delta z := \text{cell length}$$

$K(\omega) := \text{absorption profile}$

3) If I modulate my laser current I get

$$I_L = a_L V_g + b_L$$

Taking everything into account I should have

$$V_{pd}(V_g) = G(a_L V_g + b_L) e^{-K(a_L V_g + b_L) \Delta z}, \quad (P_d \propto I_{out})$$


Our aim is to get $K(\omega)$, so how to proceed?

Doppler - Broadened

1) Record

- Function generator signal (FG) in AC (zer. offset)
- SAS photodiode signal in DC and AC

2) Select a single period looking at FG signal



3) Cut both SAS signal outside the period




4) Take AC signal and correct with DC offset

$$AC \rightarrow AC - \bar{AC} + \bar{DC}$$

5) Linearize FG signal: a) Fit FG signal with linear fit

b) Replace FG values with known fit values

This is done because FG signal points are not "well behaved", especially for small voltage modulations



5) select regions far from resonance.


They should obey $V_{p2} = G(a_2 V_{f2} + b_2)$
 as $K \sim 0$ here

Fit these regions and find $G a_2$, $G b_2$.


Now we compute

$$K(a_2 V_{f2} + b_2) \Delta z = -\ln \frac{V_{p2}(V_{f2})}{G(a_2 V_{f2} + b_2)}$$

7) Since frequency ω is linear in V_{f2} we can now fit each Doppler broadened spectroscopy signal with 2 Gaussians.



8) We now have the four center positions in V_{f2} units, but we know their "true" frequencies/detunings we then fit




In this way we get $V_{f2} = a'_2 f + b'_2$, which relates with

$$f = a_2 V_{f2} + b_2 \Rightarrow V_{f2} = \frac{1}{a_2} f - \frac{b_2}{a_2}$$

$$\Rightarrow a'_2 = \frac{1}{a_2}, \quad b'_2 = -\frac{b_2}{a_2} \Rightarrow a_2 = \frac{1}{a'_2}, \quad b_2 = -\frac{b'_2}{a'_2}$$

9) We can now relate V_{f2} to f , so we plot



Remember, for Doppler Broadened spectroscopy we should have

$$\Delta f = \frac{\sqrt{\frac{k_B T}{m}}}{\lambda}$$

Doppler-Free

1) Measure F_0 , and signal with and without probe with the AC-DC trick.

You should then have

$$V_{p1}^1 = G(2v_{fg} + b_L) e^{-K_{DF}(2v_{fg} + b_L) \Delta z}$$

$$V_{p1}^2 = G(2v_{fg} + b_L) e^{-K_{DB}(2v_{fg} + b_L) \Delta z}$$

$K_{DF} := \text{doppler-free}$

$K_{DB} := \text{doppler-broadened}$

2) We can now simply get

$$(-K_{DF} + K_{DB}) \Delta z = \ln\left(\frac{V_{p1}^1}{V_{p1}^2}\right)$$

from 'theory' we should have

$$K_{DB} \approx \text{Gaussian} (1 - \beta \cdot \text{Lor})$$

$$\text{So } \ln\left(\frac{V_{p1}^1}{V_{p1}^2}\right) = \beta \cdot \text{Gaussian} \cdot \text{Lor}$$

3) So I can simply fit $\beta \cdot \text{Gauss} \cdot (\text{Lor} 1 + \text{Lor} 2 + \dots)$

with the different Lor representing the different SAS peaks

⚠ Initializing parameters for convergence is not so easy

